

## DMT Upgrade Checklist

DMT Serial # 100157

	Initials	Date
Property Identification Tag # <u>1000906</u>	JSD	5/8/25
Replace Serial # sticker	JSD	5/8/25
Affix QR code sticker	JSD	5/8/25
Install current software (43393A-B)	JSD	5/8/25
Connect printer and ensure it prints color	JSD	5/13/25
Perform Adjustment Protocol	JSD	5/13/25
Perform Calibration Protocol	JSD	5/13/25
Run MCS Solution with calibration	JSD	5/13/25
Complete Verification Testing	JSD	5/16/25
Linearity 1	JSD	5/13/25
Linearity 2	JSD	5/14/25
Linearity 3	JSD	5/16/25
Interference solutions (Acetone, MeOH, IPA)	JSD	5/13/25
Reviews completed and DMT ready for installation	JSD	5/23/25
Install DMT	JSD	6/6/25
Location: <u>Middleburg PD</u>	JSD	6/6/25
12V500 s/n: <u>MP3456</u>	JSD	↓
New log book	JSD	↓
Troubleshooting Guide	JSD	↓
Fill out DMT electronic log	JSD	6/9/25