

STATE OF VERMONT
DEPARTMENT OF PUBLIC SAFETY
VERMONT FORENSIC LABORATORY
P.O. BOX 47
WATERBURY, VERMONT 05676-0047



TEL: 802•244•8788
FAX: 802•241•5557
E-MAIL: DPS - cjs0lab@state.vt.us
<http://www.vfl.vermont.gov/>

MEMORANDUM

TO: LAW ENFORCEMENT, STATE'S ATTORNEYS AND DEFENDER GENERAL
FROM: AMANDA BOLDUC, ALCOHOL PROGRAM SUPERVISOR
SUBJECT: DATAMASTER DMT SOFTWARE UPDATE
DATE: JULY 13, 2015

A handwritten signature in blue ink, appearing to be "AB", is written over the "FROM" and "SUBJECT" lines of the memorandum.

The Vermont Forensic Laboratory's DataMaster Technical Services group is deploying a software update to all DataMaster DMT instruments. This software is being rolled out over the next few months. The analytical algorithm and all quality control parameters are not affected by this change. This update does contain minor changes to the user interface but the most significant change will be in the newly established networking capability. Once installed, each DMT that is connected to the DPS Intranet will now be able to remotely upload all tests to a computer at the Vermont Forensic Lab. Other changes to the software include:

- Fixing a glitch in the date field box on the reports menu
- Correcting spelling errors
- Changing the sample delivery display
- Clarifying question prompts
- Allowing all networked units to use centralized time
- Updating security access
- Limiting initiation of timers based on simulator temperatures
- Updating printing preferences