

## DPS - DMT

---

**From:** Dukette, Jeffery  
**Sent:** Thursday, March 03, 2016 2:07 PM  
**To:** DPS - DMT  
**Subject:** DMT Comm Log

FYI.

---

**From:** MacFarlane, Jonathan  
**Sent:** Thursday, March 03, 2016 12:06 PM  
**To:** Dukette, Jeffery  
**Subject:** RE: DMT Routine Performance Check

Jeff,

Thank you for the heads up as I saw yesterday it was out of service. I was unaware of the error and will correct it today.

Deputy MacFarlane

---

**From:** Dukette, Jeffery  
**Sent:** Thursday, March 03, 2016 8:34 AM  
**To:** MacFarlane, Jonathan  
**Cc:** Brooks, Philip  
**Subject:** DMT Routine Performance Check

Good morning,

I am just following up on the RPC performed on 2/26/2016. I left a voicemail with Chief Deputy Brooks because he was on duty at the time, but Deputy MacFarlane performed the RPC. It looks like an incorrect concentration was entered for the new simulator solution. The concentration entered into the DMT was 0.094. However, sim solution lot 15-09-100 should have a concentration of 0.099 g/210L. As a result, a solution change protocol will need to be performed before any breath samples can be collected. Please verify the concentration and lot # from the bottle/sticker and enter that data into the instrument. There is no need to change the sim solution itself, but there will still be a 30 minute waiting period as it is programmed into the software.

Thanks,

*Jeff Dukette*  
Forensic Chemist  
Vermont Forensic Laboratory  
45 State Drive  
Waterbury, VT 05671  
(802) 241-5275  
[Jeffery.Dukette@vermont.gov](mailto:Jeffery.Dukette@vermont.gov)