

Punch List
6/20/08

- 1) When the **stop time** field is left blank, the error message "Invalid Observation Time Press <Enter> to continue" appears instead of the standard message of a required field being left blank as it's supposed to. This had been corrected on a previous version of software, but is now incorrect again.
- 2) The grace period on the restart timer never resets if the subject test should time out. The DMT sat untouched for over an hour after the timer had been satisfied and the data entry screen timed out. It continued to prompt "Restart Timer? Yes/No". It should time out after fifteen minutes and return to the conventional start timer response. Please call for clarification.
- 3) The calculated range on the simulator solution change is incorrect. It's truncating there even though it's rounding in other parts of the software. This had been correct until now.
- 4) After a recalibration or a recertification, the previous Installation record disappears. It should be available in the download file, along with failing or incomplete installations.
- 5) When performing a Routine Performance Check, an Installation or any other protocol which requires the operator to input the simulator solution concentration **(except the Certification)** the solution concentration that is entered should become the default simulator solution concentration (external standard) and replace the previous lot number, concentration and range. Currently only a Simulator Solution Change protocol will allow an operator to change the simulator concentration. This is not correct. When performing an RPC or an Install, the solution that is entered is what will remain as the external standard until such time as another RPC, Simulator Solution Change or Installation occurs..
- 6) The strange blank message that occurs after an incomplete sample has not been remedied. Previously, the final Blank Test was being recorded in the field where a reason for refusal or a reason for incapable message should be. Per my fax with Scott on 6/2/08, this message was to have been removed. The DMT will always run a blank test after a purge cycle which occurs after the subject sample, however the result of this test should not be printed on a ticket. It is not printed after any other subject samples in the other testing sequences. As the software is currently, the DMT is printing the results of the final blank test under the reason for refusal message on the report. Please remove the final blank test from the report as we had requested on 6/2/08.