

ANNUAL PREVENTATIVE MAINTENANCE

DataMaster DMT: 100151
Location: NEWPORT PD
Calibration Date: 04/25/2017
Certification Date: 04/25/2017
Installation Date: 04/26/2017
APM Date: 04/10/2019
Performed by: DAVID PATLAK



Diagnostic Results

VERSIONS
DMT: 2.05
PIC: 2.08
Modem: 2.6
Questions: 2.1

TEMPERATURES

Sample Chamber = 48.8°C
Breath Tube = 42.6°C
Digital Sim = 34.0°C

SETTINGS

Lamp Voltage = 1.73 V
Cooler Voltage = 1.71 V
Bias Voltage = 80 V
Chopper Freq = 540 Hz

PUMP INFO

Flow Rate = 4.939 L/M

DETECTOR INFO

PUMP	ON	OFF
MAX(V)	0.4919	0.4961
MIN(V)	0.4902	0.4945

FILTER INFO

Filter	Value	Zero
Filter 1	0.495	true
Filter 2	0.803	true
Filter 3	1.475	true

CALIBRATION CHECK

Xq = 0.087 1.10%

APM PASSED

Simulator Vapor Test


Concentration = 0.101 g/210L
Lot # = 18-19-100
Range = 0.096 - 0.106
Simulator Vapor = 0.098 g/210L

RF Detection Test

Passed

Sample Acceptance Test

Passed



Performed by

 Date

04/10/2019

Reviewed by

 Date

4/18/19

DMT Serial Number #100151

Page 1 of 1

04/10/2019 9:19 AM

 1 of 2

Simulator Temperature Check

Simulator Serial Number:

MP3497

DMT Serial Number / Simulator Location:

100151 / Newport PD

NIST Traceable Thermometer Serial Number:

160294076

As Found

Item	Thermometer Temperature	Simulator Temperature	Difference	Within Tolerance?*
Guth 12V500 Simulator	33.92	33.99	0.07	Y

* 34.0 +/- 0.2 degree Celsius acceptable range

Was temperature adjustment performed?


Y / (N)

As Left (if applicable)

Item	Thermometer Temperature	Simulator Temperature	Difference	Within Tolerance?*
Guth 12V500 Simulator				

Notes:

Check Performed By:



Date:

4/10/19

Reviewed by:



Date:

4/18/19
4/18/19